







Critical Care • Patient Interfaces



#### FaceFit™ NIV masks

FaceFit is a comfortable single use adult hospital mask for use in non-invasive ventilation.

Two options are available; one option is for use with conventional dual limb breathing systems, and the second for bilevel and CPAP applications where an anti-asphyxiation valve and CO<sub>2</sub> leak port is required.

FaceFit masks are available in small, medium and large sizes to fit a wide range of patients. A size guide is included to help with making the right choice prior to opening the packaging.

#### Adjustable rear head pad

Allows for optimal fit for a wide range of patients

#### Lightweight

Improving patient tolerance

# Easily adjustable harness and quick release mechanism

Allowing for fast removal of the mask from the patient and optimum fit

Video available at www.intersurgical.com



## Anti-asphyxiation valve (FaceFit ported only)

Prevents rebreathing in the case of inadequate ventilation

#### CO<sub>2</sub> Leak port

(FaceFit ported only)
For effective CPAP or
bilevel ventilation

#### Low dead space

For efficient respiratory support

## 6mm stem connector with retainable cap

For supplementary oxygen or monitoring

#### Soft silicone seal

Ensuring optimum performance

### FaceFit NIV Mask - without anti-asphyxiation valve and without CO<sub>2</sub> leak port

For use with dual limb breathing systems

illustrated

Code	Description	Box Oty.
2257000	FaceFit, NIV mask, large adult	12
2256000	FaceFit, NIV mask, medium adult	12
2255000	FaceFit, NIV mask, small adult	12

#### FaceFit NIV Ported Mask - with anti-asphyxiation valve and CO2 leak port

For use with single limb breathing systems that do not have a  ${\rm CO_2}$  leak port

	Code	Description	Box Qty.
	2252000	FaceFit, NIV Ported Mask, with anti-asphyxiation valve and CO <sub>2</sub> leak port, large adult	12
	2251000	FaceFit, NIV Ported Mask, with anti-asphyxiation valve and CO <sub>2</sub> leak port, medium adult	12
	2250000	FaceFit, NIV Ported Mask, with anti-asphyxiation valve and CO <sub>2</sub> leak port, small adult	12

Make an enquiry



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.com www.intersurgical.com







The manufacturer Intersurgical Ltd is certified to ISO 9001:2008, ISO 13485:2003 and ISO 14001:2015

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

IS1.15 • Issue 6 03.18