



# Tube management programme

Help to provide a safe and tidy environment reducing the risk of disconnection and patient trauma







Critical Care • Accessories



## Tube management programme


When used in conjunction with Intersurgical breathing systems our full range of tube management options allow the clinician to position our products efficiently and comfortably around the patient area. These help to provide a safe and tidy environment reducing the risk of disconnection and patient trauma.

Intersurgical tube clips, 22F red safety cap, humidification chamber fill set spike cover, and tube clip hanger are designed to interface with the tube management clip or grip via a notch feature.


Tube clips are available to fit our Smoothbore and Flextube™ breathing systems with options available for 10mm, 15mm, 22mm and 30mm single and dual limb breathing systems.















Tube clips 10mm, 15mm, 22mm, 30mm









### **Tube support**







Our ergonomically designed tube support is a neat way of securely holding a breathing system close to the patient. Available in two sizes the tube support efficiently holds in place 10, 15, 22 and 30mm tubes with a smaller groove to grip monitoring lines and cables. The flexible neck is easily positioned to take the strain of the breathing system away from the patient.



Red safety cap or humidification chamber fill set spike cover




















### Tube management guide



























#### Make an enquiry









#### Make an enquiry



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.com www.intersurgical.com







The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2003 and ISO 14001:2015

Please think before you print
Save energy and paper.
If you must print this information
sheet please print it double sided.

IS1.40 • Issue 2 04.18