



STARMed CaStar Up

Openable hood for CPAP therapy



Critical Care ▪ Non-invasive ventilation ▪ **STARMed**



Quality, innovation and choice

Interact with us



www.intersurgical.com



STARMed CaStar Up openable CPAP hood

StarMed CaStar Up

Openable CaStar Up is a comfortable, versatile and lightweight CPAP patient interface designed to offer an alternative to traditional methods of CPAP delivery. The CaStar Up is well tolerated for extended periods and can reduce the risks associated with both endotracheal intubation and face mask use. It can be used in both semi-recumbent and supine positions providing an ideal option for CPAP therapy.

Advantages

- Allows rapid access to the patient during nursing procedures or in case of emergency
- One operator can easily remove the upper part of the hood and have full access to the patient
- Probes and catheters can remain in place, on the patient, connected to the lower part of the hood

Technical Features

- 22M/22F connectors compliant to ISO standard
- Fitted with integrated pressure gauge to check the presence of pressure within the hood
- Built in inflatable neck cushion for added patient comfort
- Comfortable and atraumatic underarm fastening straps
- Other fastening accessories are available if required
- Two new larger access ports on the collar for probes or catheters of 3.5-7.0mm Ø
- Patient access port with screw cap. Fitted with bi-directional anti-asphyxiation valve which opens automatically in the event of pressure failure
- Airtight sealing of the system is achieved by an elastic membrane making contact with the upper part of the torso, once the hood is inflated

Read more evidence on www.intersurgical.com

Features and benefits

- Patient access port**
With bi-directional anti-asphyxiation valve
- New large dimension accesses on the collar**
Ideal for single and multiple lumen probes or catheters
- Transparent**
For ease of patient monitoring
- Integrated pressure manometer**
- Easy removal**
For rapid patient access

CaStar Up openable CPAP hood with inflatable neck cushion, bi-directional anti-asphyxiation valve and patient access port

Custom packs are supplied on specific request of the customer

Code	Size	Neck circumference	Box Qty.
CA012XXL/2	Extra Extra Large	≥50cm	5
CA012XL/2	Extra Large	45-52cm	5
CA012L/2	Large	40-47cm	5
CA012M/2	Medium	34-41cm	5
CA012S/2	Small	27-34cm	5
CA012XS/2	Extra Small	17-27cm	5



www.intersurgical.com/info/starmed

Video available at www.intersurgical.com

Make an enquiry

IS1.31_ROW • Issue 4 03.16



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.com www.intersurgical.com

UK ▪ Ireland ▪ France ▪ Germany ▪ Spain ▪ Portugal ▪ Italy
Netherlands ▪ Sweden ▪ Denmark ▪ Lithuania ▪ Russia
Czech Republic ▪ Turkey ▪ South Africa ▪ China ▪ Japan
Taiwan ▪ Philippines ▪ USA ▪ Colombia ▪ Australia



Intersurgical Ltd is certified to ISO 14001, ISO 9001 and ISO 13485

Please think before you print
Save energy and paper.
If you must print this e-mail
please print it double sided